

Premier Solution for Automated
3D Confocal Raman Imaging

alpha300 *apyron*

Optimized Performance and Workflow

**NEXT
GENERATION**



apyron

with AutoBeam opto-mechanical components

The alpha300 *apyron* is the top of the line in WITec's series of Raman imaging microscopes. It has been developed to combine ease-of-use and ultimate capability by automating hardware control and offering pre-configured measurement routines. This streamlines the experimental workflow and generates reproducible results with unrivaled speed, sensitivity and resolution.



apyron [ἀπειρον] is derived from the ancient Greek word for infinity

Extensive automation

- Simplifies operation of the instrument.
- Requires less human input.
- Reduces potential sources of error.
- Enables complete remote control for use in environmental enclosures such as glove boxes.
- Enhances reproducibility.

True confocality and wavelength-optimized design

- Guarantees maximum throughput.
- Provides sharp images with resolution limited only by physical law.
- Yields spectral resolution down to 0.1 relative wavenumbers/cm (@633 nm).
- Supports FAST RAMAN IMAGING without sacrificing resolution.

What automation means

The alpha300 *apyron's* automation supports users during every step of the imaging process, making it the ideal Raman microscope for:

- Multi-user laboratories with varying requirements and levels of user experience.
- Researchers employing remote operation, such as in enclosed environments.
- Raman newcomers with advanced imaging requirements.
- Veteran Raman microscopists seeking the next performance benchmark.
- Industry facilities with recurring experimental situations and an emphasis on time-critical turnover.

Automation helps during every step of Raman imaging measurements

NEW

Setup

- TrueCal executes pre-configured calibration routines.
- Laser safety class 1 (or 1M) compliant.
- AutoBeam concept for automatic and motorized beam path alignment includes:
 - > Motorized 6-position objective turret: positions selected objective and compensates offsets
 - > Motorized illumination selector: switches between all microscopy illumination options
 - > Motorized laser coupler: delivers up to 6 laser wavelengths from UV to NIR
 - > AutoBeam output coupler: automatically optimizes signal & resolution and selects spectrometers
 - > Motorized calibration source: fast & automatic multi-point spectrometer calibration over the full spectral range
 - > Motorized Köhler illumination apertures: facilitate focusing and optimize contrast and homogeneity
 - > Motorized polarization modules: freely rotatable automated excitation polarizer and detection analyzer

Raman imaging


- TruePower determines absolute laser power with <0.1 mW accuracy.
- Instrument positions and repositions sample automatically with motorized and piezo-driven scanning stages.
- Focus stabilization actively maintains focus during entire measurements.
- Motorized objective turret includes offset-compensation for convenient precision; compatible with AFM/SNOM objectives.
- EasyLink handheld controller offers intuitive instrument operation.



Data processing, export and display

- Suite FIVE software wizard guides the user through the processing of Raman spectra, from background reduction through image generation.
- TrueMatch software helps to identify molecules by comparing measured spectra with existing databases.
- Suite FIVE facilitates the export and display of your data.

High spatial and spectral resolution simultaneously: CCl₄ in an emulsion



(A) 3D confocal Raman image of an emulsion, with a zoomed-in view of the inset. (B) shows the corresponding spectra: Green: alkane; blue: water; yellow: CCl₄ and oil. Due to the high spectral resolution of the spectroscopic system, the CCl₄ band at 460 relative wavenumbers (cm⁻¹) can be resolved at room temperature into three peaks, as shown in the zoom-in. These peaks reflect the isotopic variations of the molecules.







Image parameters: 200 x 200 x 20 pixels, 100 x 100 x 10 μm³ scan range, 6 ms integration time per spectrum, 532 nm excitation wavelength.




Raman image of a cut agate with various phases of silica and silicon oxides identified by their Raman spectra.





Raman image of a moisturizing shower gel with increased oil content. Red: oil; yellow: emulsifier; blue and cyan: aqueous phases.

The alpha300 *apyron* offers unrivaled data acquisition speed



The extremely high performance and speed of the alpha300 *apyron* microscope is vividly demonstrated with the analysis of a tungsten diselenide (WSe_2) sample.

(A) In about 2 minutes, a clear and informative $75 \times 75 \mu\text{m}^2$ Raman image of the WSe_2 flake was recorded. The flake consists of single layer areas (red), double layers (green) and multi-layers (blue). (B) shows the same measurement after smoothing. (C) A measurement of about 17 minutes produced an even sharper image. The integration time was 6 ms per pixel in both experiments. The increased signal to noise ratio was achieved by reducing the pixel size from 750 nm in (A) to 230 nm in (C). In (D) the photoluminescence image of the same flake is presented, corresponding perfectly with the Raman image. (E) White-light image of the sample.

WITec alpha300 series



**alpha300
access**

alpha300 R
Confocal Raman
Microscope

alpha300 Ri
Inverted Confocal
Raman Microscope

RISE™
Correlative
Raman-SEM Imaging

alpha300 apyron™
Automated Confocal Raman Miroscope



WE TAKE CARE
WITec uses environmentally friendly printed materials. While this policy is only a small contribution to a healthy environment, we at WITec believe that focusing on details can effect positive change in the world.

WITec Headquarters

WITec GmbH
Lise-Meitner-Str. 6
D-89081 Ulm . Germany
Phone +49 (0) 731 140700
Fax +49 (0) 731 14070200
info@witec.de
www.witec.de

WITec North America

WITec Instruments Corp.
130G Market Place Blvd.
Knoxville . TN 37922 . USA
Phone +1 865 984 4445
Fax +1 865 984 4441
info@witec-instruments.com
www.witec-instruments.com

WITec South East Asia

WITec Pte. Ltd.
25 International Business Park
#03-59A German Centre
Singapore 609916
Phone +65 9026 5667
shawn.lee@witec.biz
www.witec.de

WITec China

WITec Beijing Representative Office
Unit 1307A, Air China Plaza Tower 1
No. 36 Xiaoyun Road
Beijing, PRC., 100027
Phone +86 (0) 10 6590 0577
info.china@witec-instruments.com
www.witec.de/cn

WITec Japan

WITec K.K.
1-1-5 Furo-cho, Naka-ku,
Yokohama City, Kanagawa Pref. 231-0032
Japan
Phone +81 45 319 4277
info@witec.jp
www.witec.de/jp